



DECUS

PROGRAM LIBRARY

DECUS NO.	8-529
TITLE	OSCAR: AN OPERATING SYSTEM FOR COMPUTERIZED ANIMAL RESEARCH
AUTHOR	Dennis Kuch and John Platt
COMPANY	McMaster University Hamilton, Ontario, Canada
DATE	April 1972
SOURCE LANGUAGE	PAL III

DEC 2

LIBRARY



THE UNIVERSITY OF CHICAGO

LIBRARY

1000 S. EAST

CHICAGO, ILL.

60607

OSCAR
An Operating System
for
Computerized Animal Research

Programmer's Manual

by

Dennis O. Kuch

Development of the earliest version of OSCAR was the work of Robert Norman. Later versions were developed by John Platt. Development of various versions of OSCAR were supported by grants to John Platt from the United States Public Health Service and the National Research Council of Canada.

Purpose

This manual was written to provide the OSCAR user with

- (a) a ready reference to the facilities of the OSCAR system
- (b) the background necessary to program the user-supplied subroutines required for the conduct of specific experiments.

Prerequisites

This manual assumes that the reader has a working knowledge of the PAL-III assembly language and is familiar with the PDP-8 family of computers. Information about PAL and the PDP-8 may be obtained from the Digital Equipment Corporation's "Introduction to Programming", chapters one through six. A working knowledge may be obtained at any teletype console plugged into a PDP-8 series computer.

Table of Contents

Section	Page
1.1 Introduction	1
1.2 General functions	1
1.3 User-supplied Subroutines	2
1.4 Teletype Commands to the System	3
1.5 OSCAR Core Allocation	4
1.6 Run Buffers and the Clock Service Routine	4
2.0 OSCAR Subroutines for User Programming	5
2.1 Teletype Character Output Routines	6
2.2 Function Routines	7
2.3 Instructions for OSCAR Station Interactions	9
2.4 Arithmetic Routines	9
3.0 Writing the User-Supplied Subroutines	4
3.1 OSCAR Address Directory and Page zero Constants	4
3.2 GETVAR	5
3.3 START	8
3.4 RUN	8
3.5 OUTPUT	9

APPENDIX: A Sample User's Program

General Description

1.1 Introduction

The OSCAR system is designed to reduce the amount of programming needed for any specific class of experimental paradigms. In most cases less than four core pages of programming will be required from the user for his particular application. OSCAR oversees the servicing of the experimental stations (operant chambers, subjects' control panels, or other devices), manipulates and stores data for input or output, and provides subroutines for additional programming. OSCAR thus requires somewhat more programming than some available compiler languages but gives the user more flexibility in experimental design than these compilers. OSCAR also provides a degree of generality not available when the experimenter writes a completely new program for each new paradigm.

1.2 General functions

OSCAR's functions may be subdivided into six general categories as follows:

- 1) The Executive Monitor is the general overseer of the system and determines what gets done, when, and how. It is important to note that OSCAR controls the handling of the user's programming rather than the user's programming controlling OSCAR. That is, the Monitor is a complex interaction of subsystems which precludes the use of parts of OSCAR without engaging all of OSCAR.

- 2) The interrupt servicing routine determines system priorities in dealing with the input and output devices. First priority is given to a real-time monitor (initiated by pulses from the system clock) which ensures rapid servicing of the experimental stations. Other devices on the interrupt are the teletype reader/punch and the high speed reader/punch.

Inputs from the subject do not interrupt OSCAR but are scanned as required by the user's program.

3) Teletype routines provide operator-computer communication at all times.

4) Special input-output routines are provided for monitoring and control of experimental stations.

5) Mass data dumping may be accomplished by preprogrammed OSCAR routines.

6) OSCAR provides a set of subroutines for commonly needed functions. For example, a pseudorandom number generator, double precision arithmetic package and integer conversion routines are available for use by the user's subroutines.

1.3 User-supplied Subroutines

OSCAR requires four subroutines from the user which are called by the system during operation. These subroutines specify the particular paradigms, inputs, and outputs for classes of experiments and are described in detail in section 3. A brief description of each of the subroutines follows:

1) GETVAR: GETs VARIables (experimental parameters) from the operator.

2) START: initializes the experimental station and sets it up for the run by

(a) clearing and priming the run buffer, wherein are stored the experimental prarmeters and data for the station.

(b) turning on the output relays, which control the stimuli at the experimental station.

(c) putting a 1 in the first location of the run buffer which tells OSCAR to run the buffer.

3) RUN: is the principal and longest subroutine. It controls subject-environment interactions during the experiment proper and stores data.

4) OUTPUT: outputs data on command of the operator.

1.4 Teletype Commands to the System

The operator communicates with OSCAR through eleven numeric commands typed on the teletype: The operator types one of the following numbers on the teletype and OSCAR responds with the underlined message. The operator then types the number of the buffer to be manipulated followed by a period. If no period is given or an illegal buffer number is typed OSCAR ignores the command.

<u>Command</u>	<u>Function</u>
0:	<u>initialize</u> a station (executes GETVAR routine).
1:	<u>start</u> a station (executes START routine).
2:	<u>stop</u> a station.
3:	<u>output</u> data from a station (executes OUTPUT routine).
4:	not implemented (Mag tape dump in some versions).
5:	not implemented. (Calls library monitor in some versions).
6:	Produce ten inches of leader trailer on the high-speed punch.
7:	suppress teletype output, saving messages to the operator.
8:	examine station <u>buffer</u> . The operator types the buffer number followed by a period and then the buffer location he wants to look at. The contents of this location will be typed out by OSCAR. Both the location and its contents are in decimal.

9:	examine core location. The operator types the octal core location he wants to look at and OSCAR responds with the octal contents of that location. The operator then has the option of modifying the contents of the location. If he enters an octal number followed by a period this new number will replace the old contents. Only a period will effect modification of the contents and typing a space will skip the modification option.
::	restart the system.

1.5 OSCAR Core Allocation

OSCAR occupies 14 (decimal) pages of core, leaving 17 pages of core for the user-supplied program, the station run buffers, and the high-speed punch dump buffer. Core allocation is illustrated in Figure 1, which shows OSCAR occupying the lower nine pages and the upper five pages of core. The station run buffers, being above page zero, are said to reside in "upper core".

1.6 Run Buffers and the Clock Service Routine

Each station has associated with it a buffer storage area in upper core which contains experimental constants and variables for the particular station. On each clock pulse OSCAR transfers the first part of the buffer--the constants and variables to be used by the RUN subroutine--to page zero so that subroutines of the system can readily access it.

The first part of each run buffer is called a "run table" and contains the station number, subject identification, and experimental parameters. Two locations on page zero specify the size of the run table (TABSIZ is minus the number of locations in the run table) and the number of nonmodifiable constants (ABUMP specifies the number of locations which are not returned to the run buffer and are subsequently erased on page zero). The user must include values for these locations in his program. For example:

Figure 1

OSCAR CORE ALLOCATION

Operating System 000 - 2177
User-supplied Subroutines 2200 - XXXX

Station Run Buffers XXXX+1 - YYYY

Dump Buffer YYYY+1 - 6377
System messages 6400 - 6577
Double Precision Package 6600 - 7577
Loaders 7600 - 7777

<u>LOCATION</u>	<u>NAME</u>	<u>CONTENTS</u>
0126	TABSIZ,	-5 /Length of run table.
0127	ABUMP,	3 /Number of constants.
0130	RUNTAB,	0 /These three constants.
0131	IDNO,	0 /are required by OSCAR
0132	BOXN,	0 /
0133	NSR,	0 /Variables (or constants, then
0134	FR,	0 /variables) follow.

With each clock pulse and for each station, OSCAR goes through the following "clock service" routine which is noninterruptable:

- 1) Get the run table and bring it to page zero.
- 2) Execute the user's RUN routine if the run flag is set.
- 3) Return variables from run table to the upper core run buffer.
- 4) Get the run table for the next station and repeat.

In the upper core run buffer the first location is the station status word and indicates whether the station is to be run, stopped or ignored; and whether it has any messages to the operator. The station run flag in this location must be set by the user's START subroutine.

2. OSCAR Subroutines for User Programming

The following routines are addressed on page zero and may be accessed by an indirect jump to subroutine (JMS I XXXX) through these addresses. For example, a decimal-to-binary (DTB) conversion on a number may be affected by a JMS I 0107. More simply, A DTB instruction may be defined by DTB=JMS I 0107, reducing the complexity of coding in a user program. Such shorthand codes are usually included in the OSCAR address directory assembled with the user program.

These subroutines are not reentrant. Since the user supplied subroutines RUN and START may be executed as the result of a clock pulse which has interrupted GETVAR or OUTPUT, any routine called by GETVAR or OUTPUT can not be called by RUN or START, and vice-versa. The functions of the various system subroutines are such that the user should not be tempted to

violate this restriction.

There are four main classes of routines available in OSCAR and they will be discussed in the following order:

(1) Teletype character output

- (a) TYPET & TYPETS: character string type out
- (b) CRLF: carriage return and line feed
- (c) TAB: tabulate

(2) Function routines

- (a) RN: random number generator
- (b) PACK: stores data in dump buffer
- (c) GETI: a combination of TYPETS and DTB
- (d) CLEAR: sets run buffer variables to zero
- (e) HSPARM: reads in ASCII-coded parameters from the high-speed reader

(3) Instructions for OSCAR-station interactions

- (a) RELAY: output to station
- (b) READ: input from station

(4) Arithmetic routines

- (a) DTB: decimal-to-binary conversion
- (b) BTB: binary-to-decimal conversion
- (c) Single precision arithmetic
 - (i) SMULT: multiply
 - (ii) DIVIDE: divide
 - (iii) SQUARE: square
- (d) Double precision arithmetic

2.1 Teletype Character Output Routines

2.11 TYPET & TYPETS

TYPET is a character string output routine which types out ASCII-

coded characters on the teletype. Characters are stored starting at a location whose address must follow the TYPET subroutine call as in:

TYPET	/subroutine call
MESSAD	/address of message

The characters are stored two to a core location in six-bit (stripped of the beginning "2" or "3") ASCII code. The end of the character string is designated by a "00". The message "HELLO" which is composed of the ASCII characters 310, 305, 314, 314 and 317 would thus be stored at MESSAD as follows: MESSAD, 1005; 1414; 1700. TYPET always types a carriage return and line feed before typing the character string.

TYPETS is a subsection of the TYPET subroutine and does not require a message address when called. It uses an address established by the preceding TYPET or TYPETS call. In this manner the TYPET routine may be called to set the message address and succeeding TYPETS called to type out successive character strings. For example a "HELLO....GOODBYE" sequence with intervening programming could be coded as follows:

TYPET	/set the message address
MESSAD	
TYPETS	/type "HELLO"
:	/intervening programming
TYPETS	/type "GOODBYE"
MESSAD, 0000	/a zero sets the TYPETS pointer to MESSAD
1005;1414;1700	/ "HELLO" with 00 terminating the message
0717;1704;0231	/ "GOODBYE"
0500	

2.12 CRLF

Calling CRLF will execute a carriage return and line feed on the teletype.

2.13 TAB

Calling TAB causes the teletype to space over the number of times indicated in the location following the call as in

TAB

3

/the teletype will type 3 spaces

2.2 Function Routines

2.21 RN

The RN routine returns with a pseudo-random positive octal fraction between 0000 and 3777 in the accumulator.

2.22 GET1

This subroutine is a combination of the TYPETS and DTB subroutines, typing out a message and then inputting a number from the teletype. Use of the GET1 code saves space for this frequently occurring sequence of subroutines.

2.23 CLEAR

CLEAR sets all nonconstant locations in the run buffer to zero. It is normally used in the START routine to initialize the buffer.

2.24 PACK

This subroutine may be called by the RUN routine to store data in the high-speed punch buffer for later dumping on tape. The subroutine call requires three locations as follows:

PACK	/subroutine call
XXXX	/address of data to be stored
YYYY	/number of words to be stored

If XXXX were the location DATA and YYYY indicated that three data words were to be packed in the dump buffer, the call would be

PACK	/call
DATA	/data address
3	/3 words to be packed
:	/intervening programming
DATA,0	/data address
0	/
0	/

2.25 HSPARM

This subroutine reads in ASCII-coded parameters from the high-speed reader. It is normally called by GETVAR when a program requires a long

variable sequence of parameters. These parameters are stored sequentially beyond the run table of the appropriate run buffer. The number of parameters to be read must be in the accumulator when the subroutine is called. A dollar sign on the parameter tape designates the end of the parameter list. An end of file (EOF) message is typed if the dollar sign is encountered before enough parameters have been read.

2.3 Instructions for OSCAR Station Interactions

2.31 RELAY

This subroutine transfers the contents of the accumulator to the relay output buffer of an experimental station. How the image is transferred to different experimental stations will depend upon the particular system configuration.

RELAY is called with the desired output image in the accumulator and returns with a cleared accumulator. A calling sequence would be:

```
TAD STIM
RELAY
STIM, 4040                /simulus configuration
```

2.32 READ

Calling this subroutine transfers the contents of an input buffer from an experimental station to the accumulator. A standard subroutine returns with an image of the response inputs in the low order bits of the accumulator, although this will depend upon the particular system configuration.

2.4 Arithmetic Routines

2.41 DTB

This decimal-to-binary conversion routine is called with a cleared accumulator and returns with the converted, binary number in the accumulator. The number is taken from the teletype and must be single precision (less than 4096 in decimal). A space terminates the input number and a slash -- / --

will erase the full number if typed prior to the terminator.

This subroutine is generally used to pick up parameters at the start of the experiment in the GETVAR routine. It may also be accessed in combination with TYPETS through a GET1 call.

2.42 BTD

The binary-to-decimal conversion routine is called with the number to be outputted in the accumulator. BTD converts this number to decimal and outputs it on the teletype in a four character, integer format. Leading zeros are suppressed and spaces are typed out instead. As BTD uses the double precision interpretive package, it can not be called from the RUN routine.

Single Precision Arithmetic

Because the RUN routine can not use the double precision interpretive package, OSCAR has three single precision arithmetic routines available.

Two locations on page zero are used as registers for these single-precision subroutines. SMACL (Single-precision Multiply Accumulator, Low-order) is used for the low-order half of a number and SMACH for the high-order half. As this arithmetic is single-precision, SMACH is used only (1) for the result of a multiplication or square if the result exceeds 4095 (decimal) (2) for the dividend in a DIVIDE call if the dividend is double-precision.

2.43 SMULT

This multiply routine is entered with the multiplier in the accumulator and the multiplicand in the location following the subroutine call. The subroutine returns with the low-order part of the product in the accumulator and the high-order part in the page zero location SMACH. A typical calling sequence to multiply 1023 by 2 is:

TAD B	/Put the multiplicand after the call
DCA .+3	/
TAD A	/Get the multiplier in the accumulator
SMULT	/Call the subroutine
Ø	/Multiplicand is now here

A, 1023
B, 2

2.44 DIVIDE

To divide A by B this subroutine is entered with B in the accumulator and A in the page zero registers SMACL and SMACH. The subroutine returns with the quotient in the accumulator and the remainder in SMACH. A typical calling sequence is:

TAD A	
DCA SMACL	
TAD A + 1	/If A is single precision, SMACH must
DCA SMACH	/be cleared anyway
TAD B	
DIVIDE	/quotient will be in the accumulator
A, 1234	
5670	
B, 5231	

2.45 SQUARE

This routine is entered with the to-be-squared number in the accumulator. It returns with the double-precision square in SMACL and SMACH. The low-order part (SMACL) will be in the accumulator. A typical calling sequence to square A is:

TAD A	
SQUARE	/A squared will be in the accumulator and
A, 1234	/any overflow in SMACH

Double Precision Interpretive Package

Arithmetic Routines

The double-precision arithmetic routines are called by the user's GETVAR and OUTPUT subroutines to manipulate, store and type out parameters and data.

The double-precision package uses locations 40, 41, 42 and 43 on page

zero as an arithmetic register. These locations are labeled DPAC1, DPAC2, DPAC3, and DPAC4 respectively as they constitute the Double Precision Accumulator for the package. The low order part of a number is always stored in DPAC1 and higher order parts are stored in the three successive locations.

2.46 Special Instructions

(a) The double-precision package is entered by calling ENTER (or JMS I 7). All succeeding coding is interpreted by the double-precision interpreter and must therefore consist only of the following double-precision instructions. Enter with a cleared accumulator.

(b) Control returns to the main program when an EXIT (0000) instruction is encountered.

2.47 Arithmetic Memory Reference Instructions

(a) DPADD A

The call DPADD A causes the contents of locations A and A + 1 to be added to DPAC1 and DPAC2. The result is in DPAC1 and DPAC2.

(b) DPSUB A

DPSUB A causes the contents of locations A and A + 1 to be subtracted from DPAC1 and DPAC2. The result is in DPAC1 and DPAC2.

(c) DPMUY A

The call DPMUY A multiplies the contents of locations A and A + 1 by the contents of DPAC1 and DPAC2. The four precision answer is left in DPAC1 through DPAC4.

(d) DPDIV A

DPDIV is called with the divisor in DPAC1 and DPAC2. DPAC3 and DPAC4 should first be cleared by a DPMUY X instruction, prior to the division call. A typical calling sequence to divide the contents of B and B+1 by the contents

of A and A+1 is:

```

TAD B
DCA DPAC1
TAD B+1
DCA DPAC2
ENTER
DPMUY ONE           /Clear DPAC3 & DPAC4
DPDIV A
EXIT
B, 1234
   5670
A, 3210
   0012
ONE, 1
     0

```

The resulting quotient will be in DPAC1 and DPAC2 while the remainder will be in DPAC3 and DPAC4.

Note: Indirect addressing may be used with these memory reference instructions. However, an auto-indexing register will not be incremented by the call.

2.48 Operate Instructions

(a) DPSNA

This instruction performs a Skip on Nonzero ACcumulator where the "accumulator" referred to is the double precision accumulator; DPAC1 and DPAC2.

(b) SQROOT

Calling this routine takes the square root of the contents of DPAC1 and DPAC2 and leaves the result in DPAC1 and DPAC2.

(c) DPNEG

This instruction negates the contents of DPAC1 and DPAC2, leaving the result in DPAC1 and DPAC2.

(d) DPNOP

This is a double precision no operation and does nothing.

(e) DPIF

This instruction is followed by three addresses and performs an

arithmetic IF on the double precision accumulator DPAC1 and DPAC2. If the double precision accumulator is (1) negative, a jump to the first address is executed (2) zero, a jump to the second address or (3) greater than zero, a jump to the third address is executed. A typical calling sequence is:

DPIF	
XXXX	/where to go if DPAC1 & DPAC2 is negative
YYYY	/DPAC is zero; go to YYYY
ZZZZ	/DPAC is greater than zero; go to ZZZZ

2.49 Input/Output and Data Transfer Instructions

(a) DPOUT

This instruction performs a binary-to-decimal conversion on the double-precision accumulator and types out the result in the integer format prescribed by the call. For example DPOUT 4 outputs the number with 4 decimal places. The maximum number of digits is 7.

(b) DPGET A

The call DPGET A transfers the contents of A to DPAC1 and the contents of A + 1 to DPAC2. The contents of A and A + 1 are unchanged.

(c) DPPUT A

The call DPPUT A transfers the contents of DPAC1 and DPAC2 to locations A and A + 1. The contents of DPAC1 and DPAC2 are unchanged.

3. Writing the User-Supplied Subroutines

3.1 OSCAR Address Directory and Page Zero Constants

The user's source program (the symbolic listing in PAL-III code) is prefaced by a directory of the addresses of referenced OSCAR subroutines and registers. This is a standard listing which may be appended to all system programs. Using this directory the programmer need simply code a "TYPET" call when he wants a character string typed out and the assembler will translate this into a machine-language JMS I TYPE on the binary tape.

Another advantage of this directory, beyond parsimony in coding, is that changes in OSCAR require only a reassembling of the user's program with the new address directory and not a rewriting of the program.

The programmer must also supply page zero addresses for those locations of the run buffer which are to be referenced by the RUN routine. These locations constitute the run table. Two addresses, IDENT and BOX, are already provided by OSCAR and the user appends those constants and variables which are to be moved from the upper core run buffer to page zero during servicing of the RUN routine. Constants are listed first and followed by variables. The distinction in function between constants and variables is determined by the fact that only "variables" will be returned from page zero to the run buffer after use by the RUN routine. "Constants" therefore remain unmodified in the run buffer. How many parameters are moved from the run buffer to page zero and which parameters are treated as constants is determined by two locations on page zero. The location TABSI \bar{Z} contains the two's complement of the number of locations in the run table and determines how many parameters are moved from upper core to page zero. The location ABUMP contains the number of constants in the run table (including RUNTAB, IDENT, and BOX) and determines those parameters which will not be modified by the RUN routine. Constants must be listed first, followed by variables.

One variable frequently included in the run table is a location the contents of which points to the end of the run table in upper core. This pointer may then be used by the RUN routine to store data in the remaining locations of the run buffer.

3.2 GETVAR

The GETVAR subroutine is used to obtain experimental parameters from the operator. These should be restricted to parameters which vary among the treatments and could not be handled by program modifications at intervals

of weeks or months. Otherwise, the teletype paper expenditure becomes prohibitive and much operator time is wasted. Other parameters may be stored in the program and modified through use of the "9" command if and when necessary.

OSCAR requests four parameters in response to the operator's 0 (initialize) command:

```
INITIALIZE BUFFER # _____
IDNO _____
BOXN _____
DUMP? _____
```

The buffer number may be any number from 1 to the maximum number of buffers set up at system initialization. If the buffer is currently being run, OSCAR will mention this fact and ignore the command. Although the buffer number need not correspond to the station (or box) number --- BOXN --- it is usually less confusing to maintain this correspondence. The IDNO is an identification number for the subject to be run at this station. The station number is given in response to BOXN. The reply to DUMP? is 1 if you desire a mass data dump as specified by your RUN routine or 0 if you do not. This information tells OSCAR whether or not to accept data for dumping. The dumping routine itself (PACK) must be called by the user in his RUN subroutine as described below.

After OSCAR executes this dialog it jumps to the user's GETVAR subroutine. Included in this subroutine will be text strings requesting frequently changed parameters from the operator. This may be easily accomplished by calling:

```
TYPET
GMES
```

which will set the TYPET routine to the address of GETVAR messages GMES.

Now to output text and input parameters on the teletype one calls

GET1

which executes TYPETS and DTB, leaving the parameter in the accumulator.

These parameters may then be stored in the run table by

DCA I AUT7

where AUT7 (autoindex register 7, location 17 on page zero) has been set to point to the location after BOX in the upper core run table. OSCAR takes care of this address setting of AUT7. A complete GETVAR routine to pick up two parameters from the operator could be as simple as:

```

GETVAR, HLT
        TYPET
        GMES
        TAD M2
        DCA CTR
GLOOP,  GET1
        DCA I AUT7
        ISZ CTR
        JMP GLOOP
        JMP I GETVAR
M2,     -2
CTR,    0
GMES,   0
        4322; 0511; 1606           /"#REINF."
        7240; 0000
        2201; 2411; 1740           /"RATIO SIZE"
        2311; 3205; 7240
        0000

```

the appropriate run table would be

```

TABSIZ, -7
ABUMP,  5
RUNTAB, 0
IDENT,  0
BOX,    0
NSR,    0           /CONSTANTS
RATIO,  0           /
VAR1,   0           /VARIABLES
VAR2,   0           /

```

The GETVAR routine may be (and often is) more complicated than the above example. Complications usually involved manipulations of the parameters prior to storage in the run table. A more elaborate GETVAR routine is provided by the sample program GDT in the appendix.

3.3 START

After the buffer has been initialized, a "1" command will start the indicated buffer by calling the user's START subroutine. START should

- (a) clear the run buffer variables.
- (b) initialize all variables in the run table.
- (c) start the session by loading the appropriate stimulus configuration in the station output buffer.
- (d) set the run flag to one in the first location of the run buffer (RUNTAB).

The following coding is an example of a complete START subroutine which would accomplish these four goals:

	START, HLT	
(a)	CLEAR	/CLEAR BUFFER
	TAD NSR	/SET REINFORCEMENT COUNTER
	CIA	
(b)	DCA VAR1	
	TAD RATIO	/SET FIXED RATIO COUNTER
	CIA	
	DCA VAR2	
(c)	TAD STIM	/SET RELAYS
	RELAY	
(d)	IAC	/SET RUN FLAG
	DCA I RUNTAB	
	JMP I START	
	STIM, 2024	

The START routine in the appendix exemplifies a more complex example of this routine.

3.4 RUN

While the GETVAR, START, and OUTPUT subroutines are initiated by the operator through teletype commands, the RUN subroutine is called by OSCAR on every pulse from the system clock. If the run flag of a buffer is set to 1, OSCAR executes the RUN routine.

RUN controls all subject-environment interactions during the experiment. The user must provide for the following functions in his RUN sub-

routine.

(a) RUN must determine the current status of a station. OSCAR must know whether the station is currently in the middle of a trial, in a time out, whether reinforcement is being delivered, etc.

(b) depending upon what is happening at the station, RUN may then increment timers, modify stimuli, or check on the subject's responding.

(c) monitoring of the subject's responding then involves incrementing counters, checking timers, changing stimuli and taking other appropriate actions depending upon the contingencies.

(d) Data storage programming must be provided. At the end of a trial summed session data and/or trial data may be put away in the run buffer for later analysis by the OUTPUT subroutine.

(e) Mass data dumping may be required (see below).

(f) The RUN subroutine must stop the station by depositing 3 in RUNTAB.

If dumping of data on the high-speed paper-tape punch is desired, the subroutine PACK is used. As described previously (section 2.24) this subroutine call includes specification of the address of the data to be dumped and the number of data words to be dumped. The standard format punched on the tape is a record mark in the eighth channel of the first frame followed by successive data words split in half, six bits per tape frame. The low order half of the data word is punched first. Channel seven is always punched in all frames. OSCAR must be told to dump by replying with a "1" to the DUMP? question at station initialization or the PACK routine will be ineffective.

A complete RUN routine for the GDT program is presented in the appendix.

3.5 OUTPUT

The user's OUTPUT subroutine is accessed by OSCAR when the operator

types a "3" on the teletype and gives the desired buffer number. OUTPUT's functions may include

- (a) typing out trial-by-trial statistics if such data were stored during the run.
- (b) performing arithmetic operations on stored data and typing out summary session statistics.

As in the GETVAR routine text formatting is accomplished by the TYPET and TYPETS subroutines:

TYPET
OMES

where OMES is the first location of the data output text. The location ADRPNT on page zero is useful here as it points to the run table in upper core. By putting this address into an autoindexing register (AUT6) one may easily access the data stored in the upper core run buffer. A sample OUTPUT routine for the GDT program is provided in the appendix.

APPENDIX

A Sample User's Program:

GDT

a General Discrete Trials
program for behavioral research

A2

The GDT program will run a variety of procedures requiring response-contingent and/or time-contingent changes of stimuli (e.g., FR, VR, FT, VI). Six parameters are specified for each trial and these parameters are either read in on the high-speed reader during buffer initialization or are taken from the run buffer if they were previously stored there. The operator specifies the number of trials, the intertrial-interval stimulus (in decimal), the "hopper mode" (0 if fixed time, 1 if head-in-hopper time is used), whether or not the overwrite option is used (described below) and whether new trial parameters are to be read in (1 if yes, 0 if using the previous parameters).

The six trial parameters are

1	2	3	4	5	6
Trial time limit,	stimulus,	magazine time,	number of responses,	time interval,	ITI duration

1. The trial will be terminated, with reinforcement if such is specified by parameter #3, if the response and/or time criteria are not met by this time.
2. The trial stimulus depends upon the system configuration. It is the decimal equivalent of the buffer image to be output to the station.
3. Magazine time is the number of clock pulses during which the magazine will be energized.
4. The number of responses specifies the response criterion as in a FR or VR schedule.
5. The time interval (as in a VI or FI schedule) is specified in seconds.
6. This parameter specifies the duration of the intertrial interval in seconds. A zero indicates free operant.

Three other parameters may be specified for all stations during any

one run:

1. DROM determines what responses, if any, will reset the ITI if they occur during the ITI.
2. RMASK determines what responses will be looked at during the trial. GDT only runs single-response schedules.
3. REIN is the output buffer image for reinforcement.

The summary session statistics typed on the teletype in response to the "3" command are:

1. subject number
2. number of on-trial responses
3. session time
4. number of off-trial responses

There are two independent options available for mass data dumping:

1. If the operator replies "1" to the DUMP? interrogation during buffer initialization, the following data will be punched on the high-speed punch for each subject for each trial:

- a. subject number.
- b. stimulus and reinforcement/nonreinforcement.
- c. number of responses.
- d. time from start of trial to first response, double precision.
- e. precision.
- f. time from first response to end of trial, double precision.
- g. precision.

2. If the operator replies "1" to OVERWRITE during buffer initialization, the above trial-by-trial data will be stored over the trial parameters in the run buffer. This data can later be recorded on magnetic tape in some OSCAR configurations.


```
001      / OSCAR ADDRESS DIRECTORY
002      DO=JMS I 0
003      TYPET=DO 101
004      TYPETS=DO 102
005      CRLF=DO 105
006      TAB=DO 106
007      DTB=DO 107
008      BTB=DO 111
009      SMULT=DO 112
010      DIVIDE=DO 113
011      SQUARE=DO 114
012      RELAY=DO 115
013      READ=DO 116
014      RN=DO 117
015      GET1=DO 120
016      CLEAR=DO 121
017      PACK=DO 122
018      HSPARM=DO 123
019
020      AUT0=10
021      AUT1=11
022      AUT2=12
023      AUT3=13
024      AUT4=14
025      AUT5=15
026      AUT6=16
027      AUT7=17
028
029      ADRPNT=77
030      SMACL=75
031      SMACH=76
032      CONV=65
033
034      ENTER=DO 7
035      DPIF=1
036      DPSNA=2
037      DPOUT=100
038      DPNEG=4
039      SQROOT=3
040      DPAC1=40
041      DPAC2=41
042      DPAC3=42
043      DPAC4=43
044
045      OGET=371
046      ORUN=1510
047      ODAT=1710
048      OTAB=126
049      OBEG=2400
050
051      ACL=7701
052      CAM=7621
053      BSW=7002
```

054 SWP=7521
055 MQR=7501
056 MQL=7421
057
058 PAUSE

A6

```

001      /GDT 770
002      *ORUN
003      1510  2600      RUN
004      1511  2471      START
005      *OGET
006      0371  2400      GETVAR
007      0372  3063      OUTPUT
008      *OTAB
009      0126  7741      TABSIZ, -37
010      0127  0007      ABUMP, 7
011      0130  0000      RUNTAB, 0
012      0131  0000      IDENT, 0
013      0132  0000      BOX, 0
014      0133  0000      TRIAL, 0      /# TRIALS
015      0134  0000      ISTIM, 0      /ITISTIM
016      0135  0000      PCH, 0      /PHOTO CONTROL OF HOP TIME
017      0136  0000      OVERW, 0      /DATA OVERWRITES TRIAL VALUES
018      0137  0000      SESSR, 0
019      0140  0000
020      0141  0000      SESST, 0
021      0142  0000
022      0143  0000      ITR, 0      /INTERTRIAL RESPONSES
023      0144  0000
024      0145  0000      TLIM, 0      /TIME LIMIT (SEC)
025      0146  0000      STIM, 0      /TRIAL VALUES:
026      0147  0000      REINF, 0      /HOP TIME (.02 SEC)
027      0150  0000      RNES, 0      /RESP REQUIREMENT
028      0151  0000      TNES, 0      /TIME REQUIREMENT (SEC)
029      0152  0000      ITI, 0      /ITI (SEC)
030      0153  0000      RCNT, 0      /RESP CNT
031      0154  0000      STR, 0      /START TIME
032      0155  0000
033      0156  0000      RTIM, 0      /RUN TIME
034      0157  0000
035      0160  0000      RESP, 0      /CONTROL FLAGS:
036      0161  0000      SDIV, 0      /FREQ DIVIDER
037      0162  0000      SFLG, 0      /SECOND FLAG
038      0163  0000      TFLG, 0      /TRIAL COUNTER
039      0164  0000      IFLG, 0      /ITI FLAG
040      0165  0000      HFLG, 0      /HOPPER FLAG
041      0166  0000      DATAP, 0      /TRIAL POINTER
042
043      *ODAT
044      1710  3160      END
045
046      *OBEG
047      2400  7402      GETVAR, HLT
048      2401  4501      TYPET
049      2402  2434      GMES
050      2403  4520      GET1      /ASK TRIALS
051      2404  3231      DCA      TEM1
052      2405  1231      TAD      TEM1
053      2406  3417      DCA I      AUT7

```


054	2407	4520	GET1		/ASK ITI STIM
055	2410	3417	DCA I	AUT7	
056	2411	4520	GET1		/ASK REINF CONTROL MODE
057	2412	3417	DCA I	AUT7	
058	2413	4520	GET1		/ASK IF OVERWRITE
059	2414	3417	DCA I	AUT7	
060	2415	4520	GET1		/ASK IF NEW TRIAL VALUE
061	2416	7650	SNA CLA		
062	2417	5600	JMP I	GETVAR	/NO
063	2420	1231	TAD	TEM1	/YES
064	2421	3040	DCA	DPAC1	
065	2422	3041	DCA	DPAC2	
066	2423	4407	ENTER		
067	2424	3232	DPMUY	D6	
068	2425	0000	EXIT		
069	2426	1040	TAD	DPAC1	
070	2427	4523	HSPARM		
071	2430	5600	JMP I	GETVAR	
072	2431	0000	TEM1,	0	
073	2432	0006	D6,	6	
074	2433	0000		0	
075	2434	0000	GMES,	0	
076	2435	2422	2422		/TRIALS
077	2436	1101	1101		
078	2437	1423	1423		
079	2440	7240	7240		
080	2441	4000	4000		
081	2442	1124	1124		/ITISTIM
082	2443	1140	1140		
083	2444	2324	2324		
084	2445	1115	1115		
085	2446	7200	7200		
086	2447	1017	1017		/HOP MODE (1=PHOTO CONTROL)
087	2450	2040	2040		
088	2451	1517	1517		
089	2452	0405	0405		
090	2453	7200	7200		
091	2454	1726	1726		/OVERWRITE?
092	2455	0522	0522		
093	2456	2722	2722		
094	2457	1124	1124		
095	2460	0577	0577		
096	2461	4000	4000		
097	2462	2001	2001		/PARAMETERS?
098	2463	2201	2201		
099	2464	1505	1505		
100	2465	2405	2405		
101	2466	2223	2223		
102	2467	7740	7740		
103	2470	0000	0000		
104					
105	2471	7402	START, HLT		
106	2472	7240	STA		/CLEAR VARIABLE
107	2473	1320	TAD	STAB	

108	2474	1127	TAD	ABUMP	
109	2475	3010	DCA	AUT0	
110	2476	1126	TAD	TABSIZ	
111	2477	1127	TAD	ABUMP	
112	2500	3317	DCA	TEM2	
113	2501	3410	DCA I	AUT0	
114	2502	2317	ISZ	TEM2	
115	2503	5301	JMP	.-2	
116	2504	1126	TAD	TABSIZ	/SET TRIAL POINTER
117	2505	7041	CIA		
118	2506	1130	TAD	RUNTAB	
119	2507	3166	DCA	DATAP	
120	2510	1133	TAD	TRIAL	/SET TRIAL CNT
121	2511	7041	CIA		
122	2512	3163	DCA	TFLG	
123	2513	4321	JMS	NEW	/START NEW TRIAL
124	2514	7001	IAC		/SET RUN FLAG
125	2515	3530	DCA I	RUNTAB	
126	2516	5671	JMP I	START	
127	2517	0000		0	
128	2520	0130	TEM2,	RUNTAB	
129			STAB,		
130	2521	7402	NEW,	HLT	
131	2522	7240		STA	/GET TRIAL VALUES
132	2523	1166	TAD	DATAP	
133	2524	3010	DCA	AUT0	
134	2525	1360	TAD	M6	
135	2526	3317	DCA	TEM2	
136	2527	7240	STA		
137	2530	1357	TAD	TLIMS	
138	2531	3011	DCA	AUT1	
139	2532	1410	TAD I	AUT0	
140	2533	3411	DCA I	AUT1	
141	2534	2317	ISZ	TEM2	
142	2535	5332	JMP	.-3	
143	2536	1360	TAD	M6	/CLEAR TRIAL DATA
144	2537	3317	DCA	TEM2	
145	2540	3411	DCA I	AUT1	
146	2541	2317	ISZ	TEM2	
147	2542	5340	JMP	.-2	
148	2543	1361	TAD	M62	/SET FREQ DIVIDER
149	2544	3161	DCA	SDIV	
150	2545	1145	TAD	TLIM	/SET TIME LIMIT
151	2546	7041	CIA		
152	2547	3145	DCA	TLIM	
153	2550	1146	TAD	STIM	/LOAD RELAYS
154	2551	4515	RELAY		
155	2552	1360	TAD	M6	/UPDATE TRIAL POINTE
156	2553	7041	CIA		
157	2554	1166	TAD	DATAP	
158	2555	3166	DCA	DATAP	
159	2556	5721	JMP I	NEW	
160	2557	0145	TLIMS,	TLIM	
161	2560	7772	M6,	-6	

162	2561	7716	M62,	-62	
163					
164					
001					
002					
003			*OBEG+200		
004	2600	7402	RUN,	HLT	
005	2601	3162		DCA	SFLG /FREQUENCY DIVIDE
006	2602	2161		ISZ	SDIV /SEC?
007	2603	5207		JMP	. +4 /NO
008	2604	2162		ISZ	SFLG /YES, SET FLAG
009	2605	1356		TAD	N62 /RESET DIVIDER
010	2606	3161		DCA	SDIV
011	2607	1165		TAD	HFLG /HOPPER?
012	2610	7650		SNA CLA	
013	2611	5236		JMP	CITI /NO, CHECK FOR ITI
014	2612	1135		TAD	PCH /YES, PHOTO CONTROL?
015	2613	7650		SNA CLA	
016	2614	5222		JMP	. +6 /NO
017	2615	4516		READ	/YES, READ PHOTO CELL
018	2616	7112		CLL RTR	
019	2617	7012		RTR	
020	2620	7620		SNL CLA	
021	2621	5600		JMP I	RUN
022	2622	2165		ISZ	HFLG /HOP DONE?
023	2623	5600		JMP I	RUN /NO
024	2624	1152	ITIGO,	TAD	ITI /YES, ITI REQUEST?
025	2625	7450		SNA	
026	2626	5261		JMP	NEWT /NO, NEW TRIAL
027	2627	7041		CIA	/YES, START ITI
028	2630	3164		DCA	IFLG
029	2631	1356		TAD	N62 /SYNC DIVIDER
030	2632	3161		DCA	SDIV
031	2633	1134		TAD	ISTIM /LOAD ITI STIM
032	2634	4515		RELAY	
033	2635	5600		JMP I	RUN
034	2636	1164	CITI,	TAD	IFLG /ITI IN PROGRESS?
035	2637	7650		SNA CLA	
036	2640	5271		JMP	RCHK /NO, RUN TRIAL
037	2641	4516		READ	/RESP?
038	2642	0363		AND	DROM /SET FOR LEFT KEY OR LEVER
039	2643	7650		SNA CLA	
040	2644	5254		JMP	. +10 /NO
041	2645	1152		TAD	ITI /YES, RESET ITI
042	2646	7041		CIA	
043	2647	3164		DCA	IFLG
044	2650	2143		ISZ	ITR
045	2651	7410		SKP	
046	2652	2144		ISZ	ITR+1
047	2653	5600		JMP I	RUN
048	2654	1162		TAD	SFLG /SEC?
049	2655	7650		SNA CLA	
050	2656	5600		JMP I	RUN /NO
051	2657	2164		ISZ	IFLG /YES, ITI OVER?

A10

052	2660	5600		JMP I	RUN	/NO
053	2661	4760	NEWT,	JMS I	DMP5	/DUMP TRIAL DATA
054	2662	2163		ISZ	TFLG	/ALL TRIALS?
055	2663	5267		JMP	. +4	/NO
056	2664	1357		TAD	C3	/YES, END RUN
057	2665	3530		DCA I	RUNTAB	
058	2666	5600		JMP I	RUN	
059	2667	4761		JMS I	NEWS	/START NEW TRIAL
060	2670	5600		JMP I	RUN	
061	2671	2141	RCHK,	ISZ	SESST	/INC SESS TIME
062	2672	7410		SKP		
063	2673	2142		ISZ	SESST+1	
064	2674	1153		TAD	RCNT	/INC TRIAL TIME
065	2675	7640		SZA CLA		
066	2676	5303		JMP	. +5	
067	2677	2154		ISZ	STR	/LATENCY
068	2700	7410		SKP		
069	2701	2155		ISZ	STR+1	
070	2702	5306		JMP	. +4	
071	2703	2156		ISZ	RTIM	/RUN TIME
072	2704	7410		SKP		
073	2705	2157		ISZ	RTIM+1	
074	2706	1162		TAD	SFLG	/SECOND?
075	2707	7650		SNA CLA		
076	2710	5314		JMP	SRESP	/NO
077	2711	2145		ISZ	TLIM	/YES, TIME LIMIT?
078	2712	7410		SKP		/NO
079	2713	5346		JMP	HOP	/YES, END TRIAL
080	2714	4516	SRESP,	READ		
081	2715	0364		AND	RMASK	
082	2716	4762		JMS I	OFF5	/GET RESP OFFSETS
083	2717	7750		SPA SNA CLA		/RESP?
084	2720	5600		JMP I	RUN	/NO
085	2721	2153		ISZ	RCNT	/YES, ENOUGH RESPONSES?
086	2722	1150		TAD	RNES	
087	2723	7041		CIA		
088	2724	1153		TAD	RCNT	
089	2725	7710		SPA CLA		
090	2726	5600		JMP I	RUN	/NO
091	2727	7100		CLL		/YES, ENOUGH TIME
092	2730	1154		TAD	STR	
093	2731	1156		TAD	RTIM	
094	2732	3075		DCA	SMACL	
095	2733	7204		GLK		
096	2734	1155		TAD	STR+1	
097	2735	1157		TAD	RTIM+1	
098	2736	3076		DCA	SMACH	
099	2737	1356		TAD	N62	
100	2740	7041		CIA		
101	2741	4513		DIVIDE		
102	2742	7041		CIA		
103	2743	1151		TAD	TNES	
104	2744	7740		SMA SZA CLA		
105	2745	5600		JMP I	RUN	/NO

All

106	2746	1147	HOP,	TAD	REINF	/YES, REINF
107	2747	7450		SNA		
108	2750	5224		JMP	ITIGO	
109	2751	7041		CIA		
110	2752	3165		DCA	HFLG	
111	2753	1365		TAD	REIN	
112	2754	4515		RELAY		
113	2755	5600		JMP I	RUN	
114	2756	7716	N62,	-62		
115	2757	0003	C3,	3		
116	2760	3017	DMP5,	DMP		
117	2761	2521	NEWS,	NEW		
118	2762	3000	OFF5,	OFF		
119			/FOLLOWING CONSTANTS ARE CHAMBER SPECIFIC			
120	2763	0007	DROM,	7		/RESP CODE DURING ITI
121	2764	0007	RMASK,	7		/RESP CODE DURING TRIAL
122	2765	0003	REIN,	3		/REINF CODE
123						
001						
002			*OBEG+400			
003	3000	7402	OFF,	HLT		
004	3001	3216		DCA	TEM3	
005	3002	1160		TAD	RESP	/EXCLUSIVE OR
006	3003	0216		AND	TEM3	
007	3004	7041		CIA		
008	3005	7104		CLL	RAL	
009	3006	1160		TAD	RESP	
010	3007	1216		TAD	TEM3	
011	3010	0160		AND	RESP	/AND WITH OLD IMAGE
012	3011	7421		MQL		
013	3012	1216		TAD	TEM3	
014	3013	3160		DCA	RESP	
015	3014	7501		MQA		
016	3015	5600		JMP I	OFF	
017	3016	0000	TEM3,	0		
018						
019	3017	7402	DMP,	HLT		/SUMP TRIAL DATA
020	3020	1131		TAD	IDENT	
021	3021	3145		DCA	TLIM	
022	3022	4522		PACK		
023	3023	0145		TLIM		
024	3024	0013		13		
025	3025	1136		TAD	OVERW	/OVERWRITE?
026	3026	7650		SNA CLA		
027	3027	5250		JMP	RET	/NO
028	3030	1261		TAD	N7	
029	3031	1166		TAD	DATAP	
030	3032	3010		DCA	AUTO	
031	3033	1147		TAD	REINF	
032	3034	7640		SZA CLA		
033	3035	7001		IAC		
034	3036	1146		TAD	STIM	
035	3037	3410		DCA I	AUTO	
036	3040	1262		TAD	RCNTS	

A12

037	3041	3011		DCA	AUT1	
038	3042	1260		TAD	N5	
039	3043	3216		DCA	TEM3	
040	3044	1411		TAD I	AUT1	
041	3045	3410		DCA I	AUT0	
042	3046	2216		ISZ	TEM3	
043	3047	5244		JMP	-3	
044	3050	7100	RET,	CLL		/UPDATE RESP COUNT
045	3051	1153		TAD	RCNT	
046	3052	1137		TAD	SESSR	
047	3053	3137		DCA	SESSR	
048	3054	7204		GLK		
049	3055	1140		TAD	SESSR+1	
050	3056	3140		DCA	SESSR+1	
051	3057	5617		JMP I	DMP	
052	3060	7773	N5,	-5		
053	3061	7771	N7,	-7		
054	3062	0152	RCNTS,	RCNT-1		
055						
056	3063	7402	OUTPUT,	HLT		
057	3064	1077		TAD	ADRPNT	
058	3065	3016		DCA	AUT6	
059	3066	4501		TYPET		
060	3067	3116		OMES		
061	3070	4502		TYPETS		/IDENT
062	3071	1416		TAD I	AUT6	
063	3072	4511		BTD		
064	3073	7240		STA		
065	3074	1077		TAD	ADRPNT	
066	3075	1127		TAD	ABUMP	
067	3076	3016		DCA	AUT6	
068	3077	1315		TAD	N3	
069	3100	3314		DCA	TEM4	
070	3101	4502	LOOP,	TYPETS		
071	3102	1416		TAD I	AUT6	
072	3103	3040		DCA	DPAC1	
073	3104	1416		TAD I	AUT6	
074	3105	3041		DCA	DPAC2	
075	3106	4407		ENTER		
076	3107	0107		0107		
077	3110	0000		EXIT		
078	3111	2314		ISZ	TEM4	
079	3112	5301		JMP	LOOP	
080	3113	5663		JMP I	OUTPUT	
081	3114	0000	TEM4,	0		
082	3115	7775	N3,	-3		
083						
084	3116	0000	OMES,	0		
085	3117	2325		2325		/SUBJ
086	3120	0212		0212		
087	3121	4000		4000		
088	3122	1716		1716		/ON TRIAL RESPONSES
089	3123	4024		4024		
090	3124	2211		2211		

A13

091	3125	0114	0114
092	3126	4022	4022
093	3127	0523	0523
094	3130	2017	2017
095	3131	1623	1623
096	3132	0523	0523
097	3133	4000	4000
098	3134	2305	2305
099	3135	2323	2323
100	3136	1117	1117
101	3137	1640	1640
102	3140	2411	2411
103	3141	1505	1505
104	3142	4050	4050
105	3143	5660	5660
106	3144	6251	6251
107	3145	4000	4000
108	3146	1706	1706
109	3147	0640	0640
110	3150	2422	2422
111	3151	1101	1101
112	3152	1440	1440
113	3153	2205	2205
114	3154	2320	2320
115	3155	1716	1716
116	3156	2305	2305
117	3157	2300	2300
118	3160	0000	0

END,

/SESSION TIME(.02)

/OFF TRIAL RESPONSES

ABUMP	0127
ACL	7701
ADRFNT	0077
AUT0	0010
AUT1	0011
AUT2	0012
AUT3	0013
AUT4	0014
AUT5	0015
AUT6	0016
AUT7	0017
BOX	0132
BSW	7002
BTD	4511
CAM	7621
CITI	2636
CLEAR	4521
CONV	0065
CRLF	4505
C3	2757
DATAP	0166
DIVIDE	4513
DMP	3017
DMP5	2760
DO	4400
DFAC1	0040

A14

DPAC2	0041
DPAC3	0042
DPAC4	0043
DPIF	0001
DPNEG	0004
DPOUT	0100
DPSNA	0002
DROM	2763
DTB	4507
DG	2432
END	3160
ENTER	4407
GETVAR	2400
GET1	4520
GMES	2434
HFLG	0165
HOP	2746
HSPARM	4523
IDENT	0131
IFLG	0164
ISTIM	0134
ITI	0152
ITIGO	2624
ITR	0143
LOOP	3101
MQA	7501
MQL	7421
M6	2560
M62	2561
NEW	2521
NEWS	2761
NEWT	2661
N3	3115
N5	3060
N62	2756
N7	3061
OBEG	2400
ODAT	1710
OFF	3000
OFFS	2762
OGET	0371
OMES	3116
ORUN	1510
OTAB	0126
OUTPUT	3063
OVERW	0136
PACK	4522
PCH	0135
RCHK	2671
RCN1	0153
RCNTS	3062
READ	4516
REIN	2765
SEINF	0147

RELAY	4515
RESP	0160
RET	3050
RMASK	2764
RN	4517
RNES	0150
RTIM	0156
RUN	2600
RUNTAB	0130
SDIV	0161
SESSR	0137
SESST	0141
SFLG	0162
SMACH	0076
SMACL	0075
SMULT	4512
SQROOT	0003
SQUARE	4514
SRESP	2714
STAB	2520
START	2471
STIM	0146
STR	0154
SWP	7521
TAB	4506
TABSIZ	0126
TEM1	2431
TEM2	2517
TEM3	3016
TEM4	3114
TFLG	0163
TLIM	0145
TLIMS	2557
TNES	0151
TRIAL	0133
TYPET	4501
TYPETS	4502